MOCON® FIF OxyTraQ® Permeation system

A fast, easy to use multi-station Oxygen permeation testing system With optional purge rack

OxyTraQ The *OxyTraQ*® system...Engineered specifically for the testing of rigid packaging and containers

> Designed to allow for rapid testing of oxygen barriers and scavengers.

Features

- Combining with the Container Transfer Unit* (CTU) and optional Purge Rack, significantly increases testing throughput
- 6 stations each with its own independent sensor
- Reliable and repeatable test platform
- Intuitive user-friendly software for setup and operation
- Utilizes the innovative SmartPurge[™] technology

Robust software package

mocon

- Track real-time data From each station
- Store and export data
- Active graphing
- Monitor multiple tests







Container Transfer Units*

Purge Rack SmartPurge™

OXYGEN PERMEATION WITH SAMPLE CONDITIONII

G



Tel.: (+34) 902 450 160 Fax: (+34) 902 433 088 ermec@ermec.net

Visit us on the web: www.mocon.com

mocon

7500 Boone Ave. N Minneapolis, MN 55428 P. 763-493-6370 F. 763-493-6358 email: mocon@mocon.com



- Purging allows sample conditioning
- Allows transfer of multiple samples mounted on CTU's from pre-purge to OxyTraQ® test station:

Increase thoughput Avoid tying up sensor time

Configurable off-line Purge Rack System

Versatile Container Transfer Units* (CTUs) allow simple effect exchange between the Purge Rack and the *OxyTraQ*®



Easy mounting of various sample types and finish diameters

SmartPurge™ Technology

The innovative *SmartPurge* technology allows the *OxyTraQ*® system to operate at maximum efficiency while protecting both the samples and the oxygen sensors within the device.

- Ensures CTUs are mounted correctly and no leaks are present
- Maintains constant gas pressure and flow rate though the OxyTraO[®] system
- Extends sensor life and greatly enhances testing thoughput

| Detail | Specification |
|--------------------------|---|
| Container/Film Specifics | |
| Container materials | PET, HDPE, PP, LDPE, Monolayer and Multilayer, Barrier coated, etc. And permeable material |
| Minimum internal | • |
| Maximum external | 1.25" to 13" (32 to 330 mm)** |
| Width/Diameter | up to 4.75" (120 mm)** |
| Internal Volume | No limitation based on height and width specifications |
| Internal diameter | 0.75" to 3" (19 to 76 mm) |
| Minimum diameter | 3" (76.2 mm)** |
| Area of sample | |
| exposed to test gas | 3.875 in. ² (25 cm ²)** |
| , | ** On standard film module |
| | Specifics Container materials Minimum internal Maximum external Width/Diameter Internal Volume Internal diameter Minimum diameter Area of sample |

System Specifics

Footprint Overall dimensions of the 6 station unit (Not including computer,

gas supply and samples) 33"W x 8"D x 12"H

Operating conditions

Operating temperature 41°F to 113°F (+5°C to +45°C)

Humidity Maximum relative humidity

reaching the sensor 90% Non-condensing

MOCON commitment:

This system is another example of MOCON's 40 year commitment to innovation and quality in the design of permeation testing systems for barrier material assessment.

MOCON is a registered trademarks of MOCON and OxyTraQ and SmartPurge are trademarks of Plastic Technologies, Inc. * Container Transfer Unit - patent applied. Copyright © 2008 MOCON Inc.

MOCON reserves the right to change specifications without notice as part of our continuous program of product improvement.